

# Technology Platform Of Cell Line





http://www.quacell.com/
Quacell Biotechnology, Co., Ltd.,

SF, Building 3, No. 6, Shennong Road, Torch Development Zone, Zhongshan, China info@quacell.com
400-902-0506

QuaCell Biotechnology,Co., Ltd.

٦

Your Ultimate Biopharmaceutical Partner

# CONTENT

About Quacell	01-04
Cell line	05-06
QuaCell <sup>®</sup> CHO-K1Q Cell line	07-11
WayneLVPro™ HEK293 Cell line	12-14
WayneLVPro™ HEK293T Cell line	15
Ordering Information	16
Cell line Engineering GMP Cell Banking	17-18
Partners	



# ABOUT US



Quacell Biotechnology, Co., Ltd., headquartered in Zhongshan City, Guangdong-Hong Kong-Macao Greater Bay Area, is a rapidly growing global Biotechnology enterprise.Adhering to the requirements of ISO GMP regulations, Quacell focuses on the research and manufacture of cell medium, providing key raw materials and R&D services to the global biopharmaceutical and gene therapy companies through the innovative technologies and excellent product R & D capabilities. We promote the R & D of biopharmaceutical and biotechnology companies by providing the highest quality products and services for clients. Moreover, we develop treatments for human diseases, jointly promote the innovative development of the global biopharmaceutical industry, and improve the quality of human life.



## SUZHOU CULTURE MEDIUM PRODUCTION AND RESEARCH BASE

SUZHOU culture medium production and research base was officially put into operation in August 2022.



#### GUANGZHOU R&D CENTER

The R&D center was officially launched in February 2022.



# CELL LINE





Quacell Biotechnology, Co., Ltd., 400-902-0506 Your Ultimate Biopharmaceutical Partner www.quacell.com/

#### QuaCell<sup>®</sup> CHO-K1Q Cell line

#### The heroine of platform

- High Lac, NH4 & Osmo resistance
- Platform process, providing complete solutions to speed up R&D process
- High titer, the average IgG is more than 4.0 g/L, and the highest is 9.8 g/L
- Clear source, providing complete IND&BLA document support
- FTO, flexible global commercial license



In high Lac, NH4 & Osmo, the peak VCD of QUACELL CHO-K1Q cells was normal. In Control Condition, Cell density of QUACELL CHO-K1Q decreases slowly than CHO-K1.



#### The average titer was over 6.7g/L, and the highest titer reached 8.1g/L.



#### The stability test of TOP 10 clones showed that the titer of PDL cells fluctuated within 30%. TOP 10 clones is relatively stable and can be used for 200 L or 2000 L scale-up production.



Quacell Biotechnology, Co., Ltd., / 400-902-0506 Your Ultimate Biopharmaceutical Partner /

www.quacell.com/





### Case Study-Bispecific antibody

VCD of top10 minipools 1.40E+07 1.20E+07 1.00E+07 8.00E+06 **8** 6.00E+06 4.00E+06 2.00E+06 0.00E+00 5 6 10 11 12 13 14 4 7 8 9 3 Culture Time (Day) 







By comparing the cell pool status, growth curve and viability curve, Minipool No. 4, No. 5 and No. 7 were selected for monoclonal generalization. After a series of expansion and screening, monoclonal Fed-Batch culture was performed, and the expression level of the culture supernatant was detected using Octet RED96e Protein A. The results show that the expression levels of Top 10 clones have reached above 5.0g/L.

Quacell Biotechnology, Co., Ltd., / 400-902-0506 Your Ultimate Biopharmaceutical Partner /

www.quacell.com/

## Case Study ADC-IgG







Host cells: CHO-K1Q; Culture medium: CD02+CB7; Project Type:ADC; Expression scale:125 Flask; Incubation days:16 Days; Peak VCD: 28.2×10<sup>6</sup>cells/mL

### WayneLVPro<sup>™</sup> HEK293 Cell line

#### Lentivirus packaging tool

- High virus expression titers
- Clear source, to meet global filing requirements
- Commercial license in China, in line with global intellectual property regulations
- Suspension culture system, the highest cell density is close to 1x10<sup>7</sup>cells/mL



The results:

- 2. WayneLVPro<sup>™</sup> HEK293 Cell line is compatible with a variety of medium.

Quacell Biotechnology, Co., Ltd., / 400-902-0506 Your Ultimate Biopharmaceutical Partner / www.quacell.com/



1. The VCD of WayneLVPro<sup>™</sup> HEK293 Cell line is the highest close to 1.0×10<sup>7</sup> cells/mL.

#### **Technical Data**

WayneLVPro ™ HEK 293 Cell Line

Other HEK293 SUSPENSION Cell Line

Viral Production Cell Line

### Virus packaging expression

Comparing the titer in different volumes



The results :In the case of WayneLVPro HEK293 cell process scale-up, it can still ensure a high lentiviral titer.





WayneLVPro ™ HEK 293 Cell Line



Other HEK293 SUSPENSION Cell Line



Viral Production Cell Line



**Cell Line:** 1. WayneLVPro<sup>™</sup> HEK293 Cell Line

Medium: Wayne293® Transfection Medium

**Plasmid:** GFP-tagged plasmid

scope 48 hours after transfection.

transient efficiency, above 85%.

Transfection reagent: PEI

3. Viral Production Cell Line

2. Other HEK293 SUSPENSION Cell Line

Fluorescence expression of cells was observed under micro-

The results: WayneLVPro<sup>™</sup> HEK293 Cell Line has the highest

Three different cells were used for virus packaging under the same conditions, and the virus supernatant was harvested to infect HT1080 cells. After 48 hours, the virus titer was calculated according to the positive cell rate. WayneLVPro<sup>™</sup> HEK293 Cell Line has the highest viral titer at 1.3x10<sup>7</sup> TU/mL.

(Formula :  $TU/ml = (F \times Ci/V) \times D$ )



In the same condition, WayneLVPro™ HEK293 Cell line and VPC were used for AAV virus packaging. WayneLVPro™ HEK293 cells has highest AAV titer close to 8.5x10<sup>11</sup>VG/mL.

Quacell Biotechnology, Co., Ltd., / 400-902-0506 Your Ultimate Biopharmaceutical Partner / www.quacell.com/



#### WayneLVPro<sup>™</sup> HEK293T Cell line



WayneLVPro<sup>™</sup> HEK293T Cell line cells were batch cultured in Wayne293<sup>®</sup> and Peck VCD closed to 7.8×10<sup>6</sup> cells/mL.



WayneLVPro<sup>™</sup> HEK293T Cell Line and other HEK 293s were cultured in Wayne293<sup>®</sup>, lentivirus packaging was under the same conditions, and the WayneLVPro<sup>™</sup> HEK293T lentivirus titers were up to 5x10<sup>7</sup> TU/mL.

## Host cell platform support



## **Ordering Information**

Product	Cat.No.	Size	Storage	Application
Quacell <sup>®</sup> CHO K1Q Cell line	A13101	1.5x10 <sup>7</sup> cells/Vail	LN2	Production of protein
Wayne293™ Cell line	A23107	1.5x10 <sup>7</sup> cells/Vail	LN2	Production of protein
WayneLVPro™ HEK293 Cell line	A23109	1.5x10 <sup>7</sup> cells/Vail	LN2	Transient transfection and stable cell line development for protein
WayneLVPro™ HEK293T Cell line	A23110	1.5x10 <sup>7</sup> cells/Vail	LN2	Production of LV, AAV and other related virus vectors



#### License

- Original cell license
- Quacell license
- Global commercial application

#### Support

- Supporting series of products
- Platform technology
- IND & BLA document support



## CELL LINE ENGINEERING GMP CELL BANKING



#### PLATFORM ADVANTAGES



CHO-K1Q®Cell line High-flux instrument Octet RED 96, Improve screening rate Quacell Biotechnology, Co., Ltd., 400-902-0506 Your Ultimate Biopharmaceutical Partner www.quacell.com/







## 生物制药企业无与伦比的合作伙伴

Your Ultimate Biopharmaceutical Partner