

# QuaCell® Ready to Use Antifoam

QuaCell® Ready to Use Antifoam is designed for cell culture systems, offer quick and convenient defoaming and antifoaming effect during bioprocesses. QuaCell® Ready to Use Antifoam is made of WFI and 30% simethicone emulsion USP, and does not contain animal-derived component. QuaCell® Ready to Use Antifoam is sterilized with gamma irradiation and can be directly used for connecting to sterile production environment. QuaCell® Ready to Use Antifoam suit for clinical sample preparation and large-scale commercial production of drugs such as antibodies and fusion proteins.

## Advantage:

- Flexible customization: Bag size, pipeline and connector all are customizable.
- Convenient and Reliable: ready to use, excellent defoaming ability.

## Features

- The product is sterilized by gamma irradiation.
- This product can be directly connected to bioreactor by aseptic tubing welding.
- This product does not support filtration.
- Shake thoroughly before use.



## Order Information (standard products)

Product name	Cat.#	Size
QuaCell® Ready to Use Antifoam	P11001-02	0.2L/bag
QuaCell® Ready to Use Antifoam	P11001-05	0.5L/bag
QuaCell® Ready to Use Antifoam	P11001-25	2.5L/bag

