



## QuaCell® HEK293 CD Medium

Serum-free, animal origin free, chemically defined

- Wide adaptability, compatible with almost all HEK293 cells
- Rapid proliferation, support high-density culture, doubling time 16-24h
- Seamless connection, no acclimation process, direct subculture



Picture1.

Picture1: Comparing growth of HEK293 cells in Wayne293® and other commercial medium


## Ordering Information

Production	Cat.No.	Amount	Formulation
QuaCell® HEK293 CD Medium	A21001	1000mL	Liquid
	A22001	10L/100L/Customized	Powder




WEB


---


 <http://www.quacell.com/>

 Quacell Biotechnology, Co., Ltd.,

 5F, Building 3, No. 6, Shennong Road,

Torch Development Zone, Zhongshan, China

 [info@quacell.com](mailto:info@quacell.com)

 400-902-0506